

# Electrophoresis Mini Tankblot

## Description

High-quality blotter with safety cover for the use of 4 cassettes and 8 pads.

Transfer buffer volume: 2100 ml

Suitable with vertical electrophoresis chamber H10 Mini.

## Technical Data

|                        |                 |
|------------------------|-----------------|
| Safety lid             | Yes             |
| Cassette capacity      | 4               |
| Pad capacity           | 8               |
| Electrode material     | Platinum        |
| Transfer time          | 60 - 120 min    |
| Gel capacity           | 4 Mini gels     |
| Blotting area (W x D)  | 8 x 9.4 cm      |
| Transfer buffer volume | 2100 ml         |
| Cooling option         | ja              |
| Dimensions (W x D x H) | 21 x 15 x 25 cm |



Mini tankblot for electrophoresis



Compatible with:  
Vertical gel chamber H10 Mini