

# Polycarbonate Cryovial Storage Box 81 Places Low/High

## Description

Robust storage boxes made of polycarbonate in bright transparent colors with 81 places.

These boxes are designed for use in extreme temperatures.

Made of polycarbonate, they can be stored at -196°C to +121°C.

Fits into existing standard storage modules and provides space for all 1.5 ml to 2.0 ml microcentrifuge tubes and standard cryovials. Suitable for use in refrigerators and freezers down to -90°C minimum.

## Benefits

- Racks will withstand temperatures at -196°C - +121°C
- Vents and holes in both base and lid for draining
- Index numerals printed on the rack base (no confusion due to incorrect lid placement)
- Racks readily accommodate all standard microtubes and cryovials
- Crystal clear lid included
- Dimensions (L x W x H):  
High Box: 133 x 133 x 95 mm  
Low-Box: 133 x 133 x 54 mm



Cryovial Box Low with 81 places for 1.5 - 2.0 ml tubes



Cryovial Box High orange with 81 places, for 3-5 ml tubes

## Suitable Stainless Steel Storage Rack

Stainless steel racks with side access for freezers.

Suitable for all standard storage boxes with dimensions up to 136 x 136 x 55 mm (Low Box).

Ideal for storage of boxes between -20° C and -90° C.



Cryovial Box grid with 81 places

## Polycarbonate Cryovial Storage Box 81 Places Low/High

| Figure                                                                              | Description                                                                       | Color<br>Packaging Unit                        | Item No. | Dimensions<br>(L x W x H) |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------|----------|---------------------------|
|    | Polycarbonate Low-Box<br>81 Places<br>Lid inclusive<br><br>For 1.5 – 2.0 ml Tubes | 5 Colors mixed<br><br>5 Boxes per<br>packaging | 27-337   | 133 x 133 x 54 mm         |
|   | Polycarbonate High Box<br>81 Places<br>Lid inclusive<br><br>For 3 – 5 ml Tubes    | 5 Colors mixed<br><br>5 Boxes per<br>packaging | 27-338   | 133 x 133 x 95 mm         |
|  | Polycarbonate High Box<br>81 Places<br>Lid inclusive<br><br>For 3 – 5 ml Tubes    | Blue<br><br>5 Boxes per<br>packaging           | 26-110B  | 133 x 133 x 95 mm         |
|  | Polycarbonate High Box<br>81 Places<br>Lid inclusive<br><br>For 3 – 5 ml Tubes    | Green<br><br>5 Boxes per<br>packaging          | 26-110G  | 133 x 133 x 95 mm         |
|  | Polycarbonate High Box<br>81 Places<br>Lid inclusive<br><br>For 3 – 5 ml Tubes    | Orange<br><br>5 Boxes per<br>packaging         | 27-110O  | 133 x 133 x 95 mm         |
|  | Polycarbonate High Box<br>81 Places<br>Lid inclusive<br><br>For 3 – 5 ml Tubes    | Red<br><br>5 Boxes per<br>packaging            | 27-110R  | 133 x 133 x 95 mm         |
|  | Polycarbonate High Box<br>81 Places<br>Lid inclusive<br><br>For 3 – 5 ml Tubes    | Yellow<br><br>5 Boxes per<br>packaging         | 27-110Y  | 133 x 133 x 95 mm         |