

Polypropylene Tube Storage Module

Description

These double-sided polypropylene work racks are very flexible and can be extended as required.

The transparent lid fits onto two elements that have been plugged together.

Very dimensionally stable and robust for temperatures up to -90°C.

Specifications

- Lid included
- Durable polypropylene for temperatures down to -90°C
- Capacity for one module:
24 x 1.5/2.0 ml tubes
16 x 0.5 ml tubes
60 x 0.2 ml tubes
- Two modules plugged together can hold a 96-well PCR plate
- Dimensions with lid:
133 x 107 x 56 mm (L x W x H)



Double sided polypropylene module

Figure	Color	Description Packaging unit
	Assorted colors (yellow, red, blue, green)	Double-sided Polypropylene storage modules PU: 4 racks & 2 lids/pack