

Horizontal Stainless Steel Cooling Shelf for Standard Storage Boxes

Description

Horizontal stainless steel shelves with side access.
With practical pull-out handle.

Ideal for storing sample boxes in freezers
at -20°C and -90°C.

Suitable for use with most standard storage boxes
with maximum dimensions of 136 x 136 x 55 mm.



Horizontal stainless steel shelf,
Example illustration: 3 places horizontal / 4 places high
(order number 27-303)

Models

Order No.	Description	Capacity	Dimensions
27-303S	3 places crosswise 3 places high	9 places	415 x 140 x 168 mm
27-303	3 places crosswise 4 places high	12 places	415 x 140 x 225 mm
27-303L	3 places crosswise 5 places high	15 places	415 x 140 x 225 mm
27-304	4 places crosswise 4 places high	16 places	550 x 140 x 225 mm
27-305	4 places crosswise 5 places high	20 places	550 x 140 x 280 mm

More models on request.