

UV Handlamps - Compact Series

Description

The UV compact lamps are available in short wave (254 nm), long wave (365 nm) or in a combination of short wave and long wave.

The combined version is supplied with a switchable wavelength selector to switch between the wavelengths.

The lamps are portable and made of durable and lightweight plastic. Optional with a practical lamp stand. Also compatible with the mini darkroom C-10. All models are CE-certified.

Compatible with

- Lamp Stand J-124: Item No.: 110.4031
- Mini Viewing Cabinet C-10: Item No.: 110.5300



Technical Data

Model	Wavelegh / Tubes	Watt	Dimensions (W x D x H)	Intensity (uW/cm ³)	Part No.
UVG-11	254 nm	4 Watt	198 x 71 x 53 mm	1120	110.4025
UVGL-15*	254/365 nm 1 Split tube	4 Watt	198 x 71 x 53 mm	254 (KW) / 365 (LW)	110.4026
UVGL-25	254/365 nm, 2 tubes 1 x 254 nm / 1 x 365 nm	4 Watt	198 x 71 x 53 mm	254 (KW) / 365 (LW)	110.4028
UVL-21	365 nm	4 Watt	198 x 71 x 53 mm	1000	110.4027
UVL-23	365 nm, BLB Financial lamp	4 Watt	198 x 71 x 53 mm		110.4029

*Model UVGL-15 is not suitable for the J-124 stand due to the side slider.

Lamp stand J-124

Not suitable for model UVGL-15
Dimensions: 63 x 190 x 124 mm
Order No.: 110.4031



Mini UV viewing cabinet C-10

Order No.: 110.5300

