

Mini UV Hand Lamps

Description

Mini four-watt UV Lamps are battery operated and are excellent as a portable lamp for field use. Photo shows UVG-4, UVSL-14P and UVL-4 lamps (left to right).

Applications

- Field inspections
- Industrial inspections
- Laboratory
- Lapidary

Features

- Models include longwave, shortwave and longwave/shortwave combination
- Four Watts
- Manufactured of impact-resistant plastic which resists fading or cracking
- Lamps are small enough to fit in a pocket
- A nylon strap is attached for easy carrying
- 4 AA batteries (order separately)
- The UVSL-14P includes a slider plate for selecting either the longwave or shortwave



Mini UV lamp
UV L-4F

Specification

Model	Wavelength	Watt	Dimensions (L x T x H)	Part Number
UVG-4	254 nm	4 Watt	159 x 25 x 56 mm	110.4050
UVL-4	365 nm	4 Watt	159 x 25 x 56 mm	110.4051
UVSL-14P*	254 / 365 nm	4 Watt	159 x 25 x 56 mm	110.4053
UVL-4F	365 nm / flashlight	4 Watt	201 x 51 x 32 mm	110.4052

*Model UVSL-14P includes a slider plate for selecting either the longwave or shortwave.