

Mini UV Viewing Cabinets

Description

The mini UV darkrooms consist of a sturdy metal base with a plastic cover.

Mini Chromato-Vue Viewing Cabinets provide a compact darkroom environment for viewing materials. All cabinets feature a contrast control filter that absorbs UV energy to protect the eyes. Each cabinet is designed to be used with specific UV lamp models, which fit snugly on top of the cabinet.

Special Model

The C-15G includes a UVGL-58 multi-band 6-watt UV lamp with 254 nm and 365 nm. Please order the matching UV lamps for the models of the C-10 series separately (see models).

You can also find the matching lamps in our product range. (see page 3)

Applications

- Mineralogy
- Lab Fluorescence
- Field Investigation

Specifications

- The darkroom base panel is removable.
- Fixture for attaching a suitable lamp.
- Contrast filter absorbs UV radiation to protect the eyes.
- Large viewing window with a soft rubber inlet allows sufficient space even when wearing safety goggles.
- Light access curtain shields the booth from outside light.



C-10 Mini viewing cabinet
Basic unit without lamp



C-10E4 Mini UV viewing cabinet
with 4 Watt UV lamp EL series

Mini UV Viewing Cabinets

Models

Model Viewing Cabinet	Figure Viewing Cabinet	Compatible UV Hand Lamps Please order separately	Order No. UV Viewing Cabinet
C-10 Darkroom without hand lamp		UV Compact lamps: 4 Watt UV lamps with grip: 6 Watt	110.5300
C-10P Darkroom without hand lamp		UV hand lamp with battery Models: UVS-26P, UVSL-26P, UVL-26P	110.5302
C-10E4 Darkroom without hand lamp		EL series: 4 Watt lamps	110.5301
C-10E6 Darkroom without hand lamp		EL series: 6 Watt lamps	110.5308
C-15G Darkroom without hand lamp		Set: Darkroom C-10 6 Watt lamp UVGL-58, Shortwave/Longwave, 230 V Provides an integrated lighting system.	110.5303

UV_027_E 08/2025 Subject to technical changes and errors.

Mini UV Viewing Cabinets

Technical Data

Material base	Metal
Housing	Plastic
Viewing panel	Soft rubber edge with UV-blocking contrast control filter
Curtain	Light material, blocks ambient light
Interior dimensions (W x D x H)	229 x 254 x 178 mm
Exterior dimensions (W x D x H)	229 x 267 x 305 mm
Weight	1.6 kg (viewing cabinet without UV lamp)

Compatible Lamp Models

The lamps are available in different wavelengths.



UV Compact lamps: 4 Watt
Suitable for darkroom C-10



UV hand lamps with battery
Suitable for darkroom C-10P



UV Lampen with grip: 6 Watt
Suitable for darkroom C-10



UV-Lampen Serie EL
4 Watt, suitable for
darkroom C-10E4
6 Watt, suitable for
darkroom C-10E6