

Original Tubes for MyGo Cyclers

0.1 ml Single Tubes

Description

0.1 ml single tubes with attached lid.
Suitable for all MyGo cyclers.

Specifications:

- Low profile
- Clear
- PCR performance tested

Free of:

- Dnase/Rnase
- DNA/RNA
- Pyrogene

Packaging unit:

1000 pieces per package



MyGo single tubes with attached lid

0.1 ml 8-Strip Tubes

Description

0.1 ml Real-Time PCR tubes. The strip tubes are suitable for all MyGo-Cyclers.

Models: MyGo Pro, MyGo ESR, MyGo Mini

Specifications:

- Low profile
- Clear
- PCR performance tested

Free of:

- Dnase/Rnase
- DNA/RNA
- Pyrogene

Available in two versions:

- Strip tubes with attached cap:
Part No. 102.1504
- Strip tubes with extra cap:
Part No. tubes: 102.1502
Part No. caps 102.1503

Packaging unit:

120 strips or 120 caps per pack



Strip tubes with attached caps
(Order No. 102.1504)



8-strip tubes with extra cap
Order No. caps 102.1503,
Order No. tubes 102.1502