

0.1 ml PCR Strip Tubes for RotorGene Cycler

Description

0.1 ml Real-Time PCR Tubes.

The Strip Tubes are suitable for all RotorGene cyclers.

Models: RG-2000/3000/6000/Q

Tubes and lids are packed separately (see. fig)

Certified free of:

- RNase/DNase
- Pyrogene
- DNA/RNA

Packaging unit:

250 x 4 strips + 250 x 4 lids

Bulk, in total for 1000 reactions

Part number: 102.0170



RotorGene strip tubes



RotorGene tube lids

0.2 ml Single Tubes for RotorGene Cycler

Description

0.2 ml thin wall single tubes with attached flat cap.

Tubes are suitable for all RotorGene cyclers.

Models: RG-2000/3000/6000/Q

Certified free of:

- RNase/DNase
- Pyrogene
- DNA/RNA

Packaging unit:

1000 single tubes with attached lid

Part number: 102.0150



RotorGene single tubes