

Micro Tubes MCT 0.65 / 1.5 / 2.0 ml

Description

Microcentrifuge tubes are designed for maximum performance in a wide range of applications.

The tubes are molded from medical grade virgin polypropylene and are exceptionally clear and resistant to stretching and bulging during demanding centrifugation procedures.

The seal design allows repeated cap closing and will not pop open during boiling. All tubes feature graduation marks and writing surfaces.


Tubes are autoclavable and are supplied RNase, DNase, and Pyrogen Free.

Specifications

- Snap-Lock cap assures efficient seal, yet is still easy to 'flip' open
- Graduation marks are molded in place for optimum accuracy
- Certified Purity ensures tubes are RNase, DNase and pyrogen-free
- Ultra smooth Tube bottoms deliver perfect clarity with no cuts or scratches
- Will withstand up to 18.000 x g, in fully supported rotors



Micro Tubes MCT



Technical Data

Description	Content/Package	Packaging unit
0.65ml Micro tubes, clear, bulk, non-sterile	1000 tubes	10 pcs.
1.5ml Micro tubes, clear, bulk, non-sterile	500 tubes	10 pcs.
2.0ml Micro tubes, clear, bulk, non-sterile	500 tubes	8 pcs.