

UV Replacement Tubes VL-Series

Specifications

- Original replacement tube Vilber Lourmat
- Basic equipment quality
- for optimal results



Technical Data

Watt	Wavelengths	Part No. LTF	Model	Tubelength* ca.
6 Watt	254 nm	110.0130	T-6.C	211 mm
6 Watt	365 nm	110.0131	T-6.L	211 mm
6 Watt	White light	110.0150	T-6.WL	211 mm
8 Watt	254 nm	110.0133	T-8.C	288 mm
8 Watt	312 nm	110.0132	T-8.M	288 mm
8 Watt	365 nm	110.0134	T-8.L	288 mm
8 Watt	White light	110.0135	T-8.WL	288 mm
15 Watt	254 nm	110.0137	T-15.C	435 mm
15 Watt	312 nm	110.0136	T-15.M	435 mm
15 Watt	365 nm	110.0138	T-15.L	435 mm
15 Watt	White light	110.0153	T-15.WL	435 mm
15 Watt	Day light	110.1050	L15W/11PL-860	435 mm
20 Watt	254 nm	110.0141	T-20.C	586 mm
20 Watt	312 nm	110.0140	T-20.M	586 mm
20 Watt	365 nm	110.0139	T-20.L	586 mm
30 Watt	254 nm	110.0143	T-30.C	892 mm
30 Watt	312 nm	110.0159	T-30.M	892 mm
30 Watt	365 nm	110.0142	T-30.L	892 mm
40 Watt	254 nm	110.0145	T-40.C	1197 mm
40 Watt	312 nm	110.0146	T-40.M	1197 mm
40 Watt	365 nm	110.0144	T-40.L	1197 mm
Starter ST-151 for tubes 4 up to 15 Watt		110.0161	ST-151	-
Starter ST-111 for tubes 30 / 40 Watt		110.0162	ST-111	-

*Tubelengths without connection pins