

## 8 Strip-Tubes, Domed Cap, attached

### Description

Uniform thin-walled design for optimum heat transfer. This makes the tubes especially suitable for PCR.

Available in pre-sterile, sorted by colour, sorted by pre-sterile or transparent non-sterile.


Attached domed cap minimizes cross contamination. Each lot undergoes stringent inspection and individual lot testing ensures this tubes are certified.

#### All tubes are certified free of:

- DNA
- DNase & RNase
- Endotoxine (Pyrogene)



0.2 ml strip-tubes with domed lid, attached



### Specifications

Description	8 strip tubes 0.2 ml PCR with attached caps			
Material	Virgin polypropylene			
Compatibility	Compatible with most Thermal Cycler products on the market			
Tube Type	NP, transparent 102.5200 (3428.8.X)	NP, farblich sortiert 3428.8A.X	NP, transparent, Pre-Sterile 102.5201 (3728.8S.X)	NP, farblich sortiert Pre-Sterile (3728.8AS.X)
Quantity	Bulk, tubes & lids			
Packaging	10 bags with 12 tubes/pack ; 10 Packungen/Karton			
Autoclavable	Autoclavable at 121°C for 15 minutes at 15 PSI			
Production Standards and Processes	Product is assembled under a rigorous quality control system in a RNase/DNase free white room in accordance with ISO 9001:2000 standards. Biotix personnel are capped, masked, gloved and gowned in the assembly area. RNase/DNase free certification is verified for every lot as required.			